

Rounded Nylon bristles (SOFT) for superior cleaning
Ergonomic handle made from 100% Corn Starch (Non GMO!)

Disposal Instructions:

1. Break off toothbrush head.
2. Discard thoughtfully. Check with local recycling programs if suitable (Nylon is recyclable).
3. Place the rest of the toothbrush into your compost.

Typically takes approximately 1 year to break down into a chalk like stick in a home compost. (90 days to break down in a commercial compost.)
Congrats! You have made the right choice - any way you dispose of this, you can be sure that these toothbrushes will not be land-fill for millions of years to come - unlike many toothbrushes in use today.

Please note: The Jack N' Jill Bio Brush cannot be sterilised in boiling water or microwave. [See care & handling for details.](#)

REPLACE TOOTHBRUSH EVERY 8 WEEKS OR EARLIER